Thermofluid / *Heat Transfer Fluids*

Inspired by **temperature**

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Huber heat transfer fluids have the best thermodynamic and environmentally friendly characteristics. The choice of the correct fluid is decisive for efficient heat transfer. The working temperature range is generally the most important criterion for choosing a particular fluid. The use of the correct heat transfer fluid guarantees reliable and safe operation, optimum temperature control and a long life. The heat transfer fluids in this brochure have been chosen for use with Huber temperature control units.

If you require help in choosing the heat transfer fluid, we are happy to help:

Hotline Sales: +49 781 9603 123

Hotline Technical Support: +49 781 9603 244

*Safety data sheets and technical information for our heat transfer fluids can be found under **www.huber-online.com***

Tip for Unistats:

For Unistats, we recommend the heat transfer fluid "DW-Therm" which has a high specific heat capacity and so contributes towards an optimum heat transfer. In addition, DW-Therm covers a large temperature range, and has an impressively long service life.

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Drain tap for easy filling with thermofluid directly from the canister.

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
The table shows an indicative overview of the heat transfer fluids that may be used in each unit. When choosing the heat transfer fluid, the operating temperatures of the heat transfer fluid and unit must be considered.

DW-Therm M90.200.02
 DW-Therm HT P20.340.32
 SilOil P20.275.50
 SilOil M20.195/235.20
 SilOil M40.165/220.10
 SilOil M60.115/200.05
 SilOil M80.055/03
 SilOil M80.100/250.03
 SynOil M90.055/170.02
 MinOil M10.120.08
 MEG
 Wasser / Water

| Unistat Temperiersysteme / Temperature Control Systems | DW-Therm M90.200.02 | DW-Therm HT P20.340.32 | SilOil P20.275.50 | SilOil M20.195/235.20 | SilOil M40.165/220.10 | SilOil M60.115/200.05 | SilOil M80.055/03 | SilOil M80.100/250.03 | SynOil M90.055/170.02 | MinOil M10.120.08 | MEG | Wasser / Water |
|---|--------------------------|------------------------|-------------------|-----------------------|-----------------------|-----------------------|-------------------|-----------------------|-----------------------|-------------------|-----|----------------|
| Unistat Petite Fleur, Grande Fleur, Tango – 430w | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Unistat 510 – 540w | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Unistat 610 – 640w | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Unistat 645 – 680w | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Unistat 705 – 825w | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Unistat 904 – 950w | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Unistat 1005 – 1015w | auf Anfrage / on request | | | | | | | | | | | |
| Unistat T305 – T402, TR401 – TR402, Chili | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Unimotive 08w – 030w | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Umwälzkühler / Chillers | | | | | | | | | | | | |
| Piccolo | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Minichiller | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Unichiller 007 – 025 | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Unichiller P007 – P025 | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Unichiller 017T – 500T | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Unichiller 050 – 230, P050 – P100w | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| RotaCool | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Eintauchkühler / Immersion Cooler TC45 – TC100 | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Bad- und Umwälzthermostate / Bath Circulators | | | | | | | | | | | | |
| Einhänge-Thermostate / Immersion Circulators | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Badthermostate, Polycarbonat / Bath Circulators, Polycarbonate | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Bad-/Umwälzthermostate, Edelstahl / Bath Circulators, Stainless Steel | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Visco-Thermostate / Visco Baths | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Brücken-Thermostate / Bridge Circulators | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Kältethermostate / Cooling Circulators | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Ministat | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Variostat | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Specials | | | | | | | | | | | | |
| Bier-Forciertest-Thermostat / Beer Force-Ageing-Test-Bath | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Hotbox | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |
| Heat Transfer Station | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● | ● |

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DW-Therm is a mixture of isomeric triethoxysilanes and has been developed for hydraulically sealed systems.

We recommend the heat transfer fluid DW-Therm for Unistats. DW-Therm has a high specific heat capacity and so contributes towards an optimum heat transfer. In addition, DW-Therm covers a large temperature range, and has an impressively long service life.

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


DW-Therm HT is a mixture of partially hydrogenated naphthalenes. It is especially suited for high temperature applications using hydraulically sealed thermostats.

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


SilOil P20.275.50 is a low viscosity silicon oil, that due to its special characteristics is preferred as a cooling as well as a heating medium which can be used in temperature control units as well as in heating systems.

SilOil P20.275.50 may be used within the range 20 °C to 275 °C in connection with externally closed systems (Unistats). It can be used in open bath systems, however should be noted that P20.275.50 when used at high temperatures, in the presence of an oxidising medium such as air, or substances which act as a catalyst such as acids, alkalis, and various metallic compounds, can alter chemically. In the presence of oxidation products, an increase in viscosity or even "gelling" of the oil through a polymerisation reaction is to be expected, while in contact with catalysing substances mostly a depolymerisation process and a resulting lowering of viscosity can be observed.

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


P20.275.50 is chemically inert. Contact with concentrated alkaline solutions, oxidising acids, chlorine gas and hydrochloric acid should be avoided, particularly at elevated temperatures. Incompatible with silicone rubber.

Please pay attention to our safety data sheets and technical data sheets before use. These documents can be downloaded from www.huber-online.com or requested directly from us.

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SilOil M20.195/235.20 is a low-viscosity silicone fluid which, as a result of its special property profile, is particularly suitable for use as a cold and heat transfer medium in cryostats, thermostats and heat transfer installations.

SilOil M20.195/235.20 can be used in the range from -20 °C to 195 °C (for open systems) and to 235 °C in externally sealed systems (Unistats). It should be borne in mind that, at high temperatures, SilOil M20.195/235.20 can be chemically altered by oxidising media, such as air, or substances with a catalytic effect, such as acids, lyes and various metal compounds. An increase in viscosity, and possibly even gelling of the fluid owing to crosslinking reactions, must be expected in the presence of oxidising agents, while contact with products having a catalytic effect usually induces a process of depolymerisation, resulting in a drop in viscosity.

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


SilOil M40.165/220.10 is a low-viscosity silicone fluid which, as a result of its special property profile, is particularly suitable for use as a cold and heat transfer medium in cryostats, thermostats and heat transfer installations.

SilOil M40.165/220.10 can be used in the range from -40 °C to 165 °C and to 220 °C in externally sealed systems (Unistats). It should be borne in mind that, at high temperatures, SilOil M40.165/220.10 can be chemically altered by oxidising media, such as air, or substances with a catalytic effect, such as acids, lyes and various metal compounds. An increase in viscosity, and possibly even gelling of the fluid owing to crosslinking reactions, must be expected in the presence of oxidising agents, while contact with products having a catalytic effect usually induces a process of depolymerisation, resulting in a drop in viscosity.

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


SilOil M60.115/200.05 is a low-viscosity silicone fluid which, as a result of its special property profile, is particularly suitable for use as a cold and heat transfer medium in cryostats, thermostats and heat transfer installations.

SilOil M60.115/200.05 can be used in the range from -60 °C to 115 °C (for open systems) and to 200 °C in externally sealed systems (Unistats). It should be borne in mind that, at high temperatures, SilOil M60.115/200.05 can be chemically altered by oxidising media, such as air, or substances with a catalytic effect, such as acids, lyes and various metal compounds. An increase in viscosity, and possibly even gelling of the fluid owing to crosslinking reactions, must be expected in the presence of oxidising agents, while contact with products having a catalytic effect usually induces a process of depolymerisation, resulting in a drop in viscosity.

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


SilOil M80.055.03 is a low-viscosity silicone fluid which, as a result of its special property profile, is particularly suitable for use as a cold and heat transfer medium in cryostats, thermostats and heat transfer installations.

The fluid can be used in the range from -80 °C to 55 °C. It should be borne in mind that, at high temperatures, SilOil M80.055.03 can be chemically altered by oxidising media, such as air, or substances with a catalytic effect, such as acids, lyes and various metal compounds. An increase in viscosity, and possibly even gelling of the fluid owing to crosslinking reactions, must be expected in the presence of oxidising agents, while contact with products having a catalytic effect usually induces a process of depolymerisation, resulting in a drop in viscosity.

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


SilOil M80.100/250.03 is a low-viscosity silicone fluid which, as a result of its special property profile, is particularly suitable for use as a cold and heat transfer medium in cryostats, thermostats and heat transfer installations.

SilOil M80.100/250.03 can be used in the range from -80 °C to 100 °C (for open systems) and to 250 °C with externally sealed systems (Unistats). It should be borne in mind that, at high temperatures, SilOil M80.100/250.03 can be chemically altered by oxidising media, such as air, or substances with a catalytic effect, such as acids, lyes and various metal compounds. An increase in viscosity, and possibly even gelling of the fluid owing to crosslinking reactions, must be expected in the presence of oxidising agents, while contact with products having a catalytic effect usually induces a process of depolymerisation, resulting in a drop in viscosity.

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


Low viscous, linear polydimethylsiloxane particularly suitable as environmentally friendly cold- and heat transfer media.

SilOil M90.055/170.02 can be used in the range from -90 °C to 55 °C (for open systems) and to 170 °C in connection with externally sealed systems (for Unistats).

SilOil M90.055/170.02 is characterised by a unique combination of low viscosity and a relatively high flash point. In contrast to competitive products, SilOil M90.055/170.02 oil is not classified as dangerous good and hazardous substance. It allows an economic use even at low temperatures due to its low viscosity. Since the SilOil M90.055/170.02 is based on polydimethylsiloxane this material is easily miscible with other media on the basis of polydimethylsiloxane. SilOil M90.055/170.02 oil is colorless, odorless and non-toxic.

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


SynOil is a low viscosity Poly-alpha-olefin based fluid with a strongly enhanced low temperature characteristic. SynOil has excellent thermal stability as well as a very low tendency to evaporate. The fluid has a high viscosity index, which enables the product to flow, even at low temperatures, and not to build a disruptive film at high temperatures.

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
P20.190.40 is a high quality mineral oil based heat transfer fluid. It meets the technical safety requirements and tests according to DIN 4754 as well as the recommendations of VDI 3033 (Construction, operation and maintenance of heat transfer systems).

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


Heat transfer fluid with frost and corrosion protection characteristics for use in technical heating and cooling systems. As basis for the freezing protection, Monoethylene Glycol (MEG) is used, which due to its very high boiling point minimises loss through evaporation.

The MEG is diluted according to the table on the right side with water. Pure MEG is unsuitable for temperature control applications with Huber units.

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A protective additive for prevention and fighting the build up of algae in open bath units. It is used with water and prevents the build up of algae, bacteria and micro organisms by preventing germination.

| Bestellnummer <i>Order Number</i> | Menge <i>Quantity</i> |
|--------------------------------------|--------------------------|
| 6172 | 0,1 Ltr. |

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Corrosion is the reaction of a material with chemicals from its environment (attacking medium). The material is thereby attacked, changed, and in the worst case, destroyed. Corrosion in a watery solution is caused by an electrochemical reaction (DIN 50900-1, chapter 2.1, page 4)

Corrosion can come in various forms, depending on the cause. The forms relevant here are surface corrosion, pitting corrosion, crevice corrosion, stress corrosion/cracking and aeration corrosion.

Surface corrosion:

- Description: A slow and even change/removal of the surface.
- Cause: e.g. too small a pH-value in the water circuit.
- Effects: Optical changes of the surface, e.g. becoming dull
- Danger potential: Not dangerous, as by adding wall thickness it can be counteracted.

Pitting corrosion:

- Description: Locally increased corrosion without mechanically being stressed.
- Cause: Generally due to a higher concentration of Chloride-ions in the water circuit and the hence local destruction of the protective oxide layer on the surface.
- Outcome: Pinhole like or crater like changes in the surface.
- Danger potential: Critical as sudden leakage can occur.

Crevice corrosion:

- Description: Locally increased corrosion in narrow crevices.
- Cause: Due to inhibited diffusion, there exists a lack of Oxygen. In these areas the pH-value can thus be much reduced and metal ions from the surface become free. (Occurs for example in plate heat exchangers.) Depends on construction.
- Outcome: Changes of the material in very difficult to access areas.
- Danger potential: If at all possible, late identification of the corrosion, as the corroded area is either not visible at all or very difficult to see.

Aeration corrosion:

- Description: Locally increased corrosion under the water level.
- Cause: Differing concentrations of Oxygen in deeper layers of the water brought about by only partly filled baths and lack of circulation.
- Outcome: Local corrosive changes of the surface.

Wasser als Thermofluid und Kühlwasser / Water as a heat transfer fluid and cooling water

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Water is often used as a temperature control medium e.g. in temperature control units. In order to minimise corrosion and the thus resulting damage to the components, the water used must meet particular quality characteristics. These characteristics are described below, as well as the forms of corrosion which can be produced in areas of the unit.

Water quality

When water should be used as the temperature regulating fluid in temperature control units, then any dissolved materials, in particular ions, must be removed in order to minimise corrosion. Ions can be removed from the water in different ways. The quality of the cleaned water is then measured generally according to the electrical conductivity (in $\mu\text{S}/\text{cm}$) or respectively the electrical resistance (in $\text{M}\Omega\cdot\text{cm}$). In addition the TOC value (Total Organic Carbon) and the bacterial count can be given. There are differing standards and definitions of cleaned water, which depending on application, differ little from each other.

The following table clearly shows the various qualities of cleaned water and their classification according to various definitions.

HPW (“Highly Purified Water“) or analysis water is made using double distillation or through a continuous method such as reverse osmosis and electro dialysis.

The highest quality of water is ultra pure water. It is produced by reverse osmosis or (less often) by multiple distillations in combination with other cleaning processes, such as active carbon filters, ion exchangers, gas removal or electrical deionisation.

Quality characteristics of various cleaned waters, ASTM = American Society for Testing and Materials. CAP = College of American Pathologists

| Bezeichnung Name | Leitfähigkeit / Conductivity [$\mu\text{S}/\text{cm}$] | TOC [$\mu\text{g}/\text{mL}$] | Keimzahl / Bacterial count [KBE/mL] |
|--|---|------------------------------------|--|
| E-Wasser / De-mineralised water | < 20 | < 200 | < 1000 |
| Reinwasser ASTM Typ III, CAP Typ 3 / Pure water ASTM Type III, CAP Type 3 | < 5 | < 100 | < 100 |
| Analysewasser ASTM Typ II, CAP Typ 2 / HPW or Analysis water ASTM Type II, CAP Type 2 | < 1 | < 50 | < 10 |
| Reinstwasser ASTM Typ I, CAP Typ 1 / Ultra pure water ASTM Type I, CAP Type I | < 0,1 | < 10 | < 5 |

Wasser als Thermofluid und Kühlwasser / Water as a heat transfer fluid and cooling water

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The importance of Chloride ions in water based heat transfer fluids

Chloride ions are the main source of corrosion when using water as a heat transfer fluid in contact with metallic components. They are especially responsible for critical corrosion forms such as pitting corrosion and stress corrosion cracking. Chloride ions are adsorbed on to the (steel) surface and destroy or prevent the building of a protective (= passivating) oxide layer on the material.

Possible sources of Chloride entry:

- With E-water (or demineralised water), the exhaustion of capacity of the ion exchangers, or incorrect servicing of the facility
- Water from reverse osmosis facilities with faulty membranes, can contain Chloride ions
- The use of a biocide containing Chloride ions
- The use of highly Chlorinated sterilised drinking water, which in addition to the Chloride ions can also contain difficult to degrade organic Chlorine compounds
- In companies in which Hydrochloric Acid is used: Entry of HCl gas from the air
- Entry of salts from the environmental air, e.g. salt aerosol close to coasts

Maintaining water quality:

- If only tap water is available, the quantity of Chloride ions must be kept as low as possible in order to minimise corrosion.
- If a water circuit is operated below 50 °C and under the influence of light, then a (Chloride free!) biocide can be used to prevent algae growth.
- If there is entry of impurities or organic compounds from the air into the water circuit, then the water must be regularly changed.
- The use of corrosion inhibitors is possible, but requires careful and knowledgeable application, and is thus not recommended.
- Standard stainless steel alloys such as V2A (St 1.4301) are not resistant against corrosion caused by Chloride ions. Resistant high alloy steels such as X 8 CrNiMo 27 5 (St 1.4460) are considerably more expensive and thus mostly not economically useable.
- It is not possible to "neutralise" Chloride ions. If an acidic Chloride containing solution is for example neutralised by addition of Sodium Carbonate solution, the concentration of Chloride ions does not change, and the corrosive effects remain.

Simple checking of the water quality:

- Measurement of the pH-value; Desired range: pH = 6 – 8, ideal pH = 7 (neutral)

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- *Measurement of the electrical conductivity: desired value <math>< 50 \mu\text{S}/\text{cm}</math>*
- *Checking of colour, smell and particle content: clear/colourless, without smell or visible particles/cloudiness*


*The most economic and easiest solution when there is known Chloride content in the water is a regular water exchange with thorough cleaning using E-water (de-mineralised water)!
All comments made above are also valid for Bromide ions.*

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Corrosion resistance of stainless steel in temperature control units with differing temperatures and Chloride concentrations

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The effects of Chlorine or Chloride gases or vapours containing Chloride are to be minimised. Also Chlorine from disinfection products can cause pitting corrosion. Acidic environmental air, or water which contains a large amount of iron, or sea water also leads to corrosion. We recommend further that no untreated river or cooling tower water, which does not conform to the specifications concerning water content, be used.

Wasser als Thermofluid und Kühlwasser / Water as a heat transfer fluid and cooling water

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Use of water / Monoethylene Glycol (MEG) mixture in a temperature control unit:

The water for the mixture should have a Chloride content of < 100 mg/ltr. The specifications for other water contents should be taken from Table (page 36).

Our pure Monoethylene Glycol already contains corrosion inhibitors. The minimum mixing ratio for an adequate corrosion protection is 20 % MEG. This is equivalent to an approximate minimum working temperature of -8 °C to -10 °C. A mixing ratio of 50 % MEG should not be exceeded (viscosity).

Decalcification

With visible calcium deposits, the temperature control unit must be decalcified. For this, empty the unit and depending on the degree of calcification, make up 5 – 10 % Citric or Acetic acid solution and fill into the water circuit. Let it work for maximum 30 minutes (if possible with circulation turned on) and then empty and thoroughly rinse out all parts which had contact with the acid at least 3 times with E-water (de-mineralised water).

Dispose of in accordance with company or legal regulations. If water without any addition is used as a heat transfer fluid, then the waste water from the de-calcification can be emptied in small quantities into the drains. If the water circuit is in any way unclean, or there is a larger quantity, then dispose of it separately. In addition, take care to neutralise with Sodium Bicarbonate or Sodium Hydroxide; Check the pH-value (adjust for a pH-value of 6 – 8. Put in a collecting container for salt solutions).

Collection containers are to be clearly marked with a systematic description of their contents. The containers should be kept in a well ventilated place. These should be passed to the responsible position for waste disposal.

Taking of water samples

For the following substances, we recommend testing within 24 h after the sample was taken:

Hydrocarbonate, pH-value, Ammonia

For the following substances we recommend testing within 5 h after the sample was taken:

Free Chlorine gas, Sulphite

Our recommendation

| | |
|-------------------------------|---|
| Total hardness | 4,0 ... 8,0 °dH |
| Calcium carbonate | 0,7 mmol/L ... 1,4 mmol/L |
| pH-value: | 6,0 ... 8,5 |
| Ultra pure water, Distillate: | add 0,1 g sodium carbonate (Na ₂ CO ₃) per Liter |

Not allowed: distilled water, fully demineralised water, water with chlorid content above specified level, water which contains large amounts of iron, ammonia water, unclean water, untreated river water, sea water.

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
Corrosion resistance of brazed plate heat exchangers against water contents

The following tables give a summary of corrosion resistance of internationally used materials against various corrosives. The brazed plate heat exchangers consist of brazed stainless steel plates 1.4401 / AISI 316 as well as pure copper plates (99.9%).

There is thus the corrosion behaviour of stainless steel, Copper and the brazing material used (Copper or Nickel) to be considered. With plate heat exchanger in the standard version, Copper is used as the brazing material. The following values concerning water contents should be adhered to.

The resistance tables give a guidance value for a number of the most important chemical components. Any possibly occurring corrosion is a very complex process and caused by various contents and frequently a combination of contents together. The resistance table (table, page 36) is not complete, and serves only as a guide. The values stated are guidance values only, which under certain conditions could need to be limited further.

Temperature behaviour in relation to the corrosion resistance of plate heat exchangers with differing temperatures and Chloride concentration

The following diagram shows the allowed Chloride content for W 1.4401 / AISI 316 for plate heat exchangers with increasing temperature.

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*Plattenwärmeübertrager in Fusionstechnik komplett aus Edelstahl / *Plate heat exchangers in fusion technology completely made of stainless steel*

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Ambient Temperature Range

is the temperature range which can be attained at an ambient temperature of 20 °C by the circulator alone and with the exclusive use of electrical energy. The operating temperature, that may only be reached by using auxiliary devices, is indicated in brackets. In the case of a heating circulator the working temperature begins above room temperature (as a result of the energy introduced by the pump and the effective insulation) and ends at the upper limit of the operating temperature. The WTR of a cooling circulator begins with the lowest operating temperature of the unit and finishes with the upper temperature at which the refrigeration machine can permanently operate.

Viscosity

indicates the flow characteristics of the heat transfer fluid and it is very dependant on temperature. When operating a unit, the viscosity has a decisive influence on the temperature stability and the pump capacity.

Firepoint

is the temperature at which bath fluids continue to burn after ignition (> 5 sec). According to standard, the maximum working temperature must always be 25 °C below the fire point.

Flashpoint

is the lowest temperature at which emerging vapors can temporarily ignite (< 5 sec). The highest operating temperature specifications for bath fluids are always below the flash point and therefore exclude these risks.

Boiling Point

is the temperature at which a liquid begins to boil.

Ignition Temperature

is the lowest temperature at which the bath fluid ignites spontaneously and continues to burn without heat supply.

A series of horizontal lines for writing notes, spaced evenly down the page.

Inspired by **temperature** designed for you

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Peter Huber Kältemaschinenbau AG