

AVA & Mya 4 Control Software 2.1.2 device drivers

The table lists all the devices currently supported by AVA and Mya 4 software version 2.1.2

Device Type	Manufacturer	Model	Baud Rate	Data Bits	Parity	Stop Bits
Balance	A & D	G Series	2400	7	Even	One
	Denver	Summit	9600	8	Odd	One
	Kern	PCD 3000-2	9600	8	None	One
	Mettler	P	9600	7	Even	One
	Ohaus	Pioneer	2400	7	Even	One
	Precisa	Series XB	9600	7	Even	One
	Sartorius	Entris + M-Pact	9600	8	Odd	One
Circulator	Haake	Phoenix	9600	8	None	One
	Huber	Pilot	9600	8	None	One
	Julabo	F	4800	7	Even	One
	Lauda	Integral XT	9600	8	None	One
	Lauda	Proline	9600	8	None	One
Conductivity Meter	Knick	703	4800	7	Even	One
Gas Controller	Bronkhorst	EI-Flow	38400	8	None	One
Gas Monitor	M+W Instruments (+Bronkhorst)	Mass-Stream	38400	8	None	One
Hotplate Stirrer	Heidolph	HEI-End (+ Hei-Connect)	9600	7	Even	One
	Heidolph	Hei-Plate Mix'n'Heat Ultimate	9600	7	Even	One
	IKA	Ret Control	9600	7	Even	One
	Radleys	Carousel Advanced (+Carousel Connect)	9600	7	Even	One
pH Meter	Jenway	3520 pH	9600	8	None	One
	Knick	Portamess 913 pH	4800	7	Even	One
	Knick	Portavo 904 pH	N/A	N/A	N/A	N/A
Pressure Probe	Omega	PXM409-BGUSBH	115200	8	None	One
Peristaltic Pumps	Cole-Palmer	Masterflex	4800	7	None	One
	Heidolph	Hei-FLOW + pumpdrive PD	9600	8	None	One
	Heidolph	Hei-FLOW Precision 01	9600	8	None	One
	Heidolph	Hei-FLOW Precision 06	9600	8	None	One
	Heidolph	Hei-FLOW Ultimate 120	9600	8	None	One
	Heidolph	Hei-FLOW Ultimate 600	9600	8	None	One
	ISMATEC	ISM321 REGLO-CPF	9600	8	None	One
	KNF	SIMDOS 02	9600	8	None	One
	KNF	SIMDOS 10	9600	8	None	One
	Watson-Marlow	323Du	9600	8	None	Two
	Watson-Marlow	520Du	9600	8	None	Two
	Watson-Marlow	530Du	9600	8	None	Two
Piston Pumps	KNAUER	AZURA P 2.1S	9600	8	None	One
	KNAUER	AZURA P 4.1S	9600	8	None	One
	Teledyne ISCO	ReaXus M1	9600	8	None	One
Stirrer	Heidolph	Hei-TORQUE + RZR	9600	8	None	One
	IKA	EUROSTAR power control-visc	9600	7	Even	One
	Lara	Stirrer	9600	8	None	One
	Radleys	RS	9600	8	None	One
Syringe	Harvard	Model '11' Elite	9600	8	None	Two
	Harvard	Model '11' Plus	9600	8	None	Two
	New Era (+WPI)	NE-1000 (+ AL1000-220)	19200	8	None	One
Temp Probe	Radleys	Data Hub Temperature Probe	N/A	N/A	N/A	N/A
Vac Controller	Vacuubrand	CVC 3000	19200	8	None	One