

SHAKING



Hei-SHAKE Reciprocating Core

THE MOST STURDY WITH
36 MM STROKE



Energy Efficiency*



Safety



Limitless possibilities

Stroke of 36 mm & broad rpm range
(1 – 300 rpm)

Made for various vessels

Flexible attachments (i.e. changeable angle)
for tool-free assembly Platforms 1000
and 2000

Intuitive handling

Standby button in view and one knob with
LED ring to set values and identify the active
functions

Display and timer for an easy control of the
process parameters

RS232 for remote control and data recording
via software

Innovative body design for a secure handling

High robustness and efficiency

Excellent load to footprint ratio (max.
15 kg load)

Up to 50% lower energy consumption when
compared to other similar shakers

Low temperature rise on platform



MADE IN GERMANY

3-year warranty on all devices and an
operational lifespan of more than 10 years

For more information visit [our website](#)

*Claims based on internal tests. Not made by third party.


RECOMMENDED ACCESSORIES



Set for upgrade
to platform 2000



Separatory
funnel attachment



Attch. for test tubes & Falcons
(centrifuge tubes) (45-90°)



Multi-tier
attachment