

SHAKING



Hei-SHAKE Orbital Core

THE ALL-ROUNDER UP TO 15 KG



Energy Efficiency*



Safety



Limitless possibilities

Broad rpm range (1 – 530 rpm)

Made for various vessels

Flexible attachments (i.e. changeable angle) for tool-free assembly Platforms 1000 and 2000

Intuitive handling

Standby button in view and one knob with LED ring to set values and identify the active functions

Display and timer for an easy control of the process parameters

RS232 for remote control and data recording via software

Innovative body design for a secure handling

High robustness and efficiency

Excellent load to footprint ratio (max. 15 kg load)

Up to 50% lower energy consumption when compared to similar shakers

Low temperature rise on platform

For more information visit [our website](#)


*Claims based on internal tests. Not made by third party.



MADE IN GERMANY

3-year warranty on all devices and an operational lifespan of more than 10 years


RECOMMENDED ACCESSORIES



Set for upgrade to platform 2000



Separatory funnel attachment



Attachment for test tubes, Falcon tubes, and Eppis



Multi-tier attachment