

SHAKING



Hei-SHAKE Salto Core

THE MOST ROBUST FOR
HEAVY APPLICATIONS

DIN 19529



Energy Efficiency*



Safety



The Powerhouse of Shaking

Broad rpm range (1 – 32 rpm)

Made for various type of vessels up to 300 mm height and 136 mm diameter

Capacity for up to 12 bottles

Intuitive handling

Standby button in view and one knob with LED ring to set values and identify the active functions

Display and timer for an easy control of the process parameters

RS232 for remote control and data recording via software

High efficiency and safety

Excellent load to footprint ratio
(max. 30 kg load)

Lower energy consumption

Overload detector for immediate shut down



MADE IN GERMANY

3-year warranty on all devices and an operational lifespan of more than 10 years

For more information visit [our website](#)

*Claims based on internal tests. Not made by third party.

RECOMMENDED ACCESSORIES



Spare tension plate for caps



Attachment set