

SHAKING



Hei-SHAKE Wave Core

THE BIOSPECIALIST FOR
CONTINUOUS RINSING



Energy Efficiency*



Safety



Smooth and efficient shaking

Broad rpm range (1 – 80 rpm)

Available in 10° angle to avoid foaming

Available in 5° angle to avoid foaming,
breakage and spillage

Intuitive handling

Standby button in view and one knob with
LED ring to set values and identify the active
functions

Display and timer for an easy control of the
process parameters

RS232 for remote control and data recording
via software

High robustness and safety

Excellent load to footprint ratio
(max. 15 kg load)

Up to 50% lower energy consumption
when compared to other similar shakers

Innovative body design for a
secure handling

Low temperature rise on platform



MADE IN GERMANY

3-year warranty on all devices and an
operational lifespan of more than 10 years

For more information visit [our website](#)

*Claims based on internal tests. Not made by third party.

RECOMMENDED ACCESSORIES



Rubber mat 1000 & 2000



Set for upgrade to platform 2000



Multi-tier attachment